

Ask for economical repair and service offers!

Maintenance & Calibration Pricelist Sweden, Denmark, Norway and Iceland	Level of calibration ¹⁾			
	Price	10/1	10/2	10/3
Mechanical Single-channel	€	46	60	71
Mechanical Multichannel	€	67	86	107
Mechanical Multichannel 0,5-10 ul	€	107	125	144
Electronic Single-channel	€	60	71	83
Electronic Multichannel	€	69	88	109
Electronic Multichannel 0,5-10 ul	€	109	127	146
Mechanical/Electronic Dispenser	€	63	72	87

Extras Pricelist	Price €
FINAS Accredited certificate (ISO 17025:2005) Single-channel ²	39
FINAS Accredited certificate (ISO 17025:2005) Multichannel ²	62
As Found Performance Test Single-Channel	24
As Found Performance Test Multichannel	40
As Found Performance Test Multichannel <10µl	79
Pipette Decontamination	35
Small Parts Supplement	7
Cost Estimation	30
Delivery	35

1) Performance test measurements/volume. For example 10/3 means 10 measurements of three volumes (10%, 50% and 100% of the pipette's maximum volume).

2) An accredited calibration certificate (per ISO 17025:2005) is provided if ordered.

During a Sartorius service, testing is carried out according to ISO 8655 specifications, unless otherwise agreed with the customer. Sartorius standard tips are used in Performance testing.

Service Procedure:

- The pipette is opened and carefully cleaned.
- When necessary, worn out parts are replaced.
- Seals and o-rings are lubricated.
- Performance testing (ISO 8655:2002 specifications)

The cost of spares are added separately to the invoice. Every channel in the multichannel pipette is tested individually.

A service and performance test reports are dispatched with the serviced pipette.

Sartorius standard tips are used in performance testing. If the customer uses other brands of tips, they need to be sent in with the pipette.



Pipette performance test is carried out according to ISO 8655:2002 specifications. Sartorius Biohit Liquid Handling Oy's calibration laboratory is a FINAS accredited calibration laboratory no. K041. The Finnish Accreditation Service (FINAS) co-operates with international accreditation organisations, which co-operate to harmonise and promote the consistency of accreditation. FINAS is a member of the European co-operation for Accreditation (EA) and is part of the Multilateral Recognition Agreement (MLA, MRA).

Sartorius Biohit Liquid Handling Oy
Laippatie 1
FI-00880 Helsinki
Finland
Tel.: +358 50 4141 628

Quotation requests and orders: sales.nordic@sartorius.com
Information requests: service.nordic@sartorius.com